

Glycan Profiling of Crude Sample

Comparison of CHO and Lec1 mutant cell lines

As is well known, cells have specific glycomes depending on the origins: species, organs, differentiation stages, and/or health conditions. Therefore, glycan are usually called "the face of the cell". Up until the present, there have not been any simple and high-sensitive ways to identify cells focusing on the difference of glycomes.

In this technical note, we compare glycomes between CHO and Lec1 mutant cells, from which GlcNAc-transferase I was knocked out, using GlycoStation™. The differences in glycomes are detected from fluorescence patterns due to interactions of Cy3 labeled cell membrane glycoproteins with lectins immobilized on LecChip™.

◆ GlcNAc-transferase I

As is shown in Fig.1, GlcNAc-transferase I (GlcNAc-T I) plays important role in adding GlcNAc onto the α 1-3 and α 1-6 Man of complex and hybrid type N-glycan core structure. Therefore, in a case of GlcNAc-T I deficient Lec1 mutant, we can theoretically expect that Lec1 lacks complex and hybrid type N-glycans, resulting in increase of high-mannose type N-glycans instead.




Fig.1 GlcNAc-T I works to change glycan structures from high-mannose to complex and hybrid type N-glycans.

◆ Characteristic differences in lectin signals

We observed characteristic differences in N-glycan binding lectins from the lectin microarray analyses. The signals of complex type N-glycan binders (PHA(L), PHA(E), ACG), that for α 2,3-Sia binder (MAL I), and that for asialo binders (RCA120) decreased in Lec1 (see Fig.2(A)), whereas the signals of high-mannose binders (GNA, HHL, PWM, Calsepa, PSA, LCA) increased in Lec1 (see Fig.2(B)). These findings are quite reasonable taking into consideration that Lec1 is a GlcNAc-T I deficient mutant.




Fig. 2 Comparison of Net Cy3 fluorescence intensities between CHO and Lec1 mutant cells

(A) Lectin signals higher in CHO (B) Lectin signals higher in Lec1

◆ Dependence of signal intensity on sample concentration




Fig.3 Concentration dependent signal intensity curves for ACG and Calsepa

Quantitative differential comparison is possible in GlycoStation™, utilizing a fact that the lectin microarrays signals get stronger with increasing sample concentration. Fig.3 shows concentration dependent signal intensity curves for ACG and Calsepa using membrane fractions taken from CHO and Lec1 mutant cells. For your information, ACG typically binds to Sia α 2-3Gal β 1-4GlcNAc, and Calsepa has specificity to monnose type N-glycans. These findings promise you that the result of differential analysis doesn't change with sample concentration and the profiling pattern is also kept constant, as far as the experiment is done in a proper concentration range without signal saturation. This is one of the great appealing points of GlycoStation™.

- Protocol -

1. Sample preparation and fractionation

- 1-1. Wash cells (5×10^6) by PBA several times and store at -80°C .
- 1-2. Fractionate the cells into membrane and cytoplasm fractions by using a commercial Kit¹⁾.

2. Cy3 labeling

- 2-1. Measure protein concentration with a commercial Micro BCA Protein Assay Reagent Kit²⁾ (reaction time = 2h).
- 2-2. Prepare 20 μL sample volume with a concentration of 50 $\mu\text{g}/\text{ml}$ using PBS³⁾, then mix it with Cy3 Mono-Reactive dye 100 μg labeling⁴⁾.
- 2-3. React it for 1h in a dark place at R.T.
- 2-4. Add 300 μL TBS⁶⁾ into a gel filtration column⁵⁾, then centrifuge at 1,500 $\times g$, for 1 min at 4°C .
- 2-5. Repeat it two times.
- 2-6. Add the whole sample prepared at 2-3 and 25 μL TBS into the gel filtration column prepared at 2-5, then centrifuge at 1,500 $\times g$, for 2min, at 4°C , and remove excess Cy3.

3. Measurement

- 3-1. Dilute samples with a Probing Solution⁷⁾ (in a range from 2 $\mu\text{g}/\text{mL}$ to 31.25 ng/mL).
- 3-2. Wash LecChip™⁸⁾ by a Probing Solution three times, then apply samples into wells (100 $\mu\text{L}/\text{well}$).
- 3-3. Incubate LecChips at 20°C over night.
- 3-4. Measure fluorescence patterns without any washing of LecChips by GlycoStation™ Reader 1200⁹⁾.
- 3-5. Analyse the results by Array-Pro™ Analyzer¹⁰⁾ and GlycoStation™ Tools¹¹⁾.

Note

- 1) ProteoExtract® Subcellular Proteome Extraction Kit (Calbiochem, #539790)
- 2) Micro BCA Protein Assay Reagent Kit (PIERCE, #23235)
- 3) PBS(-) pH7.3
- 4) Use the equivalent amount of Cy3 Mono-Reactive dye pack (GE, #PA23011) to that used for 100 μg protein labeling.
- 5) Zeba Desalt Spin Columns, 0.5ml (Thermo, #89883)
- 6) TBS pH7.5
- 7) Probing Solution (GlycoTechnica)
- 8) LecChip™ (GlycoTechnica)
- 9) GlycoStation™ Reader 1200 (GlycoTechnica)
- 10) Array-Pro™ Analyzer ver.4.5 (MEDIA CYBERNETICS)
- 11) GlycoStation™ Tools (GlycoTechnica)